

INTUITIVE

Da Vinci SureForm

Achieve more with your stapler



Achieve more with your stapler

Eliminate guesswork with the intelligence and control of a SureForm stapler.¹

At its core is SmartFire technology, which monitors tissue compression before and during firing, making automatic adjustments to optimize the staple line.

A 120° cone of articulation¹ allows you the freedom to approach tissue where it naturally lies, delivering greater maneuverability around target tissue and potentially reducing tissue tension.²

When you choose SureForm staplers, you can potentially reduce uncertainty and increase stapling performance.¹

¹ Results from internal testing on file.

² Testing performed in benchtop surgeon evaluations. Surgeons using SureForm staplers applied significantly less target tissue tension during position and clamping than Echelon Flex, Endo GIA, and Signia. Data on file as of March 2019.



Complete surgeon control

All SureForm staplers are seamlessly integrated with our da Vinci surgical systems¹ to allow you full control of positioning, clamping, and firing from the surgeon console.

Intuitive offers a full portfolio of SureForm staplers to ensure you have the right instrument to meet your procedural needs.

SureForm 60



SureForm 45



SureForm 45 Curved-Tip



8 mm SureForm 30 Curved-Tip



¹ Compatible with the da Vinci X and da Vinci Xi surgical systems.

Greater range of motion

A 120° cone of articulation¹ gives you the freedom to maneuver your stapler and position it where you prefer.

01
60° Down


02
60° Left

03
60° Right


04
60° Up



Utilize a greater range of motion than the human hand or leading laparoscopic staplers with SureForm staplers.²



Powered Echelon Flex



Signia / Endo GIA Ultra

¹ Results from internal testing on file.

² Compared to Endo GIA Ultra with Tri-Staple technology, Signia with Tri-Staple 2.0 technology, and Powered Echelon Plus with GST as of March 2019.

Intelligent stapling

With one press of a pedal you activate SmartFire technology, which measures tissue compression before and during firing.

SmartFire technology uses more than 1,000 measurements per second¹ to make automatic adjustments to the firing process as staples are being formed and as the transection is being made.

This advanced firing technology is easy to use and helps ensure a more consistent staple line, all while preventing damage across a range of tissue thicknesses.²




¹ Results from internal testing on file.

² Compared to Endo GIA Ultra with Tri-Staple technology, Signia with Tri-Staple 2.0 technology, and Powered Echelon Plus with GST as of March 2019.

Optimally formed staples across a range of tissue thicknesses¹


Percentage of optimally formed staples in thick tissue (>4 mm)¹



¹ Testing performed with black reloads on in vivo porcine stomach tissue with both 45 mm and 60 mm devices. SureForm 60 has significantly more optimally formed staples than Endo GIA and Signia in thick tissue. SureForm 45 has significantly more optimally formed staples than all the three staplers in thick tissue. Data on file as of March 2019.

Bring tissue together with less outer layer tearing¹

Successful tissue approximation rate¹



¹ Testing performed with black reloads on in vivo porcine stomach tissue with both 45 mm and 60 mm devices. SureForm 45 and SureForm 60 have significantly less outer tissue tearing than Echelon Flex. Data as of March 2019.

More consistency

SmartFire's automatic measurements and adjustments during firing help you maintain control of the tissue.

This protects the integrity of the cut line and helps achieve the expected length of transection. These consistent results can help you avoid using an unanticipated reload, potentially reducing cost.

Deliver a more consistent cut line length across thin and thick tissue.¹

Figure A: Average cutline length (60 mm staplers)





Figure B: Average cutline length (45 mm staplers)



¹ Testing performed with black reloads on ex vivo porcine stomach tissue. The results of cut line consistency were determined via linear regression of the effect of tissue thickness on cut line length. The cutline length of SureForm 60 and Signia with 60 mm reloads were NOT significantly affected by tissue thickness, while the cutline length of Echelon Flex and Endo GIA with 60 mm reloads were affected by tissue thickness (Figure A). The cutline length of SureForm 45 was NOT significantly affected by tissue thickness, while the cutline length of all three other staplers with 45 mm reloads were affected by tissue thickness (Figure B). Data on file as of March 2019.


Less tension

With the full cone of articulation and anvil centric design of SureForm staplers, you can position the device exactly where you want, minimizing tension on the target tissue and stress to the surrounding anatomy.


Less tension on target tissue and less stress on surrounding tissue during positioning and clamping¹

01
Tension on target tissue.

02
Stress on surrounding tissue.




Average tension applied to target tissue while positioning and clamping stapler



Measured in newtons as a unit of force

Average stress applied to surrounding tissue while positioning and clamping stapler




Measured in megapascals as a unit of stress



¹ Testing performed in benchtop surgeon evaluations. Surgeons using SureForm staplers applied significantly less target tissue tension during positioning and clamping than Echelon Flex, Endo GIA, and Signia. Surgeons using SureForm staplers applied significantly less surrounding tissue stress during positioning and clamping than Echelon Flex, Endo GIA, and Signia. Data on file as of March 2019.

More stability during firing¹



Note: Average amount of tip movement during firing in a side-by-side benchtop test.

Images are for illustration purposes only and are not to the exact size.

¹As tested in ex vivo surgeon evaluations. SureForm 60 and competitive 60mm stapler data displayed. SureForm 45 and SureForm 60 demonstrated significantly less tip movement during firing than Echelon Flex, Signia, and Endo GIA. Data on file as of March 2019.



8 mm SureForm 30 Curved-Tip

Small stapler. Big impact.

Experience the performance of a SureForm stapler, now on a smaller scale.¹





¹ Compared to 12 mm SureForm staplers.

More maneuverability¹

01
Fits through 8 mm cannula

02
120° articulation

03
Narrower curved-tip anvil¹



The 8 mm SureForm 30 Curved-Tip fits through the standard cannula used with other da Vinci instruments, enabling you to choose the optimal angle to access tissue through any available port.

With SureForm's standard 120° of articulation, a 30 mm reload length, and a narrower¹, curved-tip anvil, this stapler arms you with a greater range of motion² in tight spaces and the flexibility to position and maneuver it as needed.

¹ Compared to 12 mm SureForm staplers.

² Compared to Echelon Powered Vascular Stapler and Signia Small Diameter Reloads.

More precision

A narrower anvil fits through smaller dissection windows¹ allowing for easy access around vessels and the 120° of articulation helps you precisely position the stapler with less tissue tension.²

Less tension on target tissue and less stress on surrounding tissue during positioning, clamping, and firing.³


36%

Less tension on target tissues²
with 8 mm SureForm 30


48%

Less stress on surrounding tissues²
with 8 mm SureForm 30

Average tension applied to target tissue¹



Average stress applied to surrounding tissue¹



■ 8 mm SureForm 30 ■ Echelon Powered Vascular Stapler ■ Signia Small Diameter Reloads

¹ Compared to 12 mm SureForm staplers.

² Compared to Ethicon Echelon Powered Vascular Stapler and Signia Small Diameter Reloads.

³ Testing performed in benchtop surgeon evaluations. Surgeons using 8 mm SureForm 30 applied significantly less target tissue tension during positioning, clamping, and firing than Echelon PVS and Signia. Surgeons using 8 mm SureForm 30 applied significantly less surrounding tissue stress during positioning, clamping, and firing than Echelon PVS and Signia. Data on file as of March 2022.

Performance you can depend on

Built on the same trusted foundation as 12 mm SureForm staplers, the 8 mm SureForm 30 is powered by SmartFire technology, which monitors tissue compression before and during firing, making automatic adjustments to optimize the staple line.

30%

Less tip motion during clamping and firing with 8 mm SureForm 30 compared to the Echelon Powered Vascular stapler

46%

Less tip motion during clamping and firing with 8 mm SureForm 30 compared to the Signia Small Diameter Reloads

More stability during firing¹




1.4 mm

**8 mm SureForm
30 Curved-Tip
tip motion**



2.0 mm

**Echelon Powered
Vascular Stapler
tip motion**



2.6 mm

**Signia Small
Diameter
Reloads
tip motion**

Average amount of tip movement during firing in a side-by-side benchtop test. Images are for illustration purposes only and are not to the exact size.

¹ As tested in ex vivo surgeon evaluations. 8 mm SureForm 30 demonstrated significantly less tip movement during firing than Echelon Powered Vascular Stapler and Signia Small Diameter Reloads. Data on file as of March 2022.

Secure hemostasis

Substantially equivalent hemostasis¹

8 mm SureForm 30 (2 row) has been tested with grey and white reloads on in vivo porcine vessels and pedicles, and demonstrated substantially equivalent hemostasis to SureForm 45 (3 row).¹

Significantly better hemostasis²

8 mm SureForm 30 has been tested with white reloads on in vivo porcine vessels and bundles, and demonstrated significantly better hemostasis than Echelon Powered Vascular Stapler and Signia Small Diameter Reloads.²

¹ Compared to SureForm 45. Testing performed with grey and white reloads on in vivo porcine vessels. 8 mm SureForm 30 demonstrated equivalent hemostasis performance to SureForm 45 at $p < 0.05$. Data on file as of April 2021.

² Compared to Echelon Powered Vascular Stapler and Signia Small Diameter Reloads. Testing performed with white reloads on in vivo porcine vessels. 8 mm SureForm 30 demonstrated significantly better hemostasis than Echelon PVS and Signia Small Diameter Reloads at $p < 0.05$. Data on file as of March 2022.

SureForm specifications

Stapler

- 8 mm SureForm 30:
32 mm staple line/30 mm cut line
- 12 mm SureForm 45:
45 mm staple line/40 mm cut line
- 12 mm SureForm 60:
60 mm staple line/55 mm cut line
- An anvil-centric design
- 12 mm SureForm 45 and SureForm 60 operate through a 12 mm cannula
- 8 mm SureForm 30 operates through an 8 mm cannula
- 12 fires per device
- Single patient use

SureForm 60 reload

White, blue, green, black reloads

Six rows of staples

New knife with every reload



SureForm 45 reload

Grey, white, blue, green, black reloads

Six rows of staples


New knife with every reload



8 mm SureForm 30 reload

Grey, white, blue reloads





Four rows of staples



SureForm staplers

Ordering information





Instruments for da Vinci X & Xi surgical systems

	Product description	Qty/ box	Uses	Part number
	12 mm SureForm 60 stapler	6	Single	480460
	12 mm SureForm 45 stapler	6	Single	480445
	12 mm SureForm 45 curved-tip stapler	6	Single	480545
	8 mm SureForm 30 curved-tip stapler	6	Single	488530

SureForm 60






Ordering information

Reloads for da Vinci X & Xi SureForm 60 staplers

	Product description	Qty/ box	Uses	Part number
	12 mm SureForm 60 white reload (2.5 mm, 6-row)	12	1 firing	48360W
	12 mm SureForm 60 blue reload (3.5 mm, 6-row)	12	1 firing	48360B
	12 mm SureForm 60 green reload (4.3 mm, 6-row)	12	1 firing	48360G
	12 mm SureForm 60 black reload (4.6 mm, 6-row)	12	1 firing	48360T

Ordering information




Reloads for da Vinci X & Xi SureForm 45 staplers

	Product description	Qty/ box	Uses	Part number
	12 mm SureForm 45 grey reload (2.0 mm, 6-row)	12	1 firing	48345M
	12 mm SureForm 45 white reload (2.5 mm, 6-row)	12	1 firing	48345W
	12 mm SureForm 45 blue reload (3.5 mm, 6-row)	12	1 firing	48345B
	12 mm SureForm 45 green reload (4.3 mm, 6-row)	12	1 firing	48345G
	12 mm SureForm 45 black reload (4.6 mm, 6-row)	12	1 firing	48345T

8 mm SureForm 30

Ordering information



Reloads for da Vinci X & Xi SureForm 30 stapler

	Product description	Qty/ box	Uses	Part number
	8 mm SureForm 30 curved-tip stapler grey reload (2.0 mm, 4-row)	12	1 firing	48230M
	8 mm SureForm 30 curved-tip stapler white reload (2.5 mm, 4-row)	12	1 firing	48230W
	8 mm SureForm 30 curved-tip stapler blue reload (3.5 mm, 4-row)	12	1 firing	48230B



SureForm 60 and 45 stapler accessories

Ordering information

Cannulas for da Vinci X & Xi SureForm 60 and 45 staplers





	Product description	Qty/ box	Uses	Part number
	12 mm and stapler cannula (100 mm)	1	NA	470375
	12 mm and stapler cannula, long (150 mm)	1	NA	470389

Disposables for da Vinci X & Xi SureForm staplers

	Product description	Qty/ box	Uses	Part number
	12 mm and stapler cannula seal	10	1 procedure	470380
	12-8 mm reducer	6	1 procedure	470381

Ordering information

Obturator for da Vinci X & Xi SureForm 60 and 45 staplers


	Product description	Qty/ box	Uses	Part number
	12 mm and stapler blunt obturator	1	NA	470376
	12 mm and stapler blunt obturator, long	1	NA	470390
	12 mm and stapler bladeless obturator	1	NA	470395
	12 mm and stapler bladeless obturator, long	1	NA	470396

**8 mm SureForm 30
stapler accessories**

**Ordering
information**

Cannulas for da Vinci X & Xi SureForm 30 stapler

	Product description	Qty/ box	Uses	Part number
--	----------------------------	---------------------	-------------	------------------------




8 mm cannula

1

NA

470002



8 mm cannula, long


1

NA

470004

Obturator for da Vinci X & Xi SureForm 30 stapler

	Product description	Qty/ box	Uses	Part number
--	----------------------------	---------------------	-------------	------------------------




8 mm blunt obturator

1

NA

470008



8 mm blunt obturator, long

1

NA

470009

8 mm SureForm 30 stapler accessories

Ordering information

Disposables for da Vinci X & Xi SureForm 30 stapler

Product description	Qty/box	Uses	Part number
12-8 mm reducer	6	1 procedure	470381¹
8 mm bladeless obturator (optical)	6	1 procedure	470359
8 mm bladeless obturator, long (optical)	6	1 procedure	470360
5-8 mm cannula seal	10	1 procedure	470361
Universal cannula seal (5-12 mm)	10	1 procedure	470500

To learn more about da Vinci SureForm stapler products, visit:

www.intuitive.com

To contact a representative or receive additional information:

Visit www.intuitive.com or call Intuitive Customer Service: +41.21.821.2020 or +800.0821.2020 (toll free)

European Regional Offices

France

Intuitive Surgical SAS
Cité de la Photonique,
bâtiment Gienah
11 avenue de Canteranne
33600 Pessac, FR
Tel: +33 5 86 07 01 01
N° de TVA : FR 64 417 889 177
SIREN 417 889 177
Capital Social 37000€
RCS Bordeaux Code APE 4646Z

United Kingdom

Intuitive Surgical Ltd
The Schrodinger Building
Oxford Science Park
Heatley Rd, Oxford, OX4 4GE, UK
Tel: +44 (0) 1865 411561
N° de TVA : GB943 2213 49
N° d'immatriculation
6727762

European Headquarters (Switzerland)

Intuitive Surgical Sàrl
1, chemin des Mûriers
1170 Aubonne, Suisse
Tél. : +41 21 821 2000
Fax : +41 21 821 2001
N° de TVA : CHE-113.407.162 VAT

Germany

Intuitive Surgical Deutschland GmbH
Am Flughafen 6
79108 Freiburg, DE
MWST #: DE200038727
Amtsgericht Freiburg, HRB 720585
Geschäftsführer:
David Chenux, Tim Heisterkamp,
Jamie Samath, Gary Loeb

¹ Version 12 or higher.

Da Vinci X & Xi Surgical Systems

The Intuitive Surgical Endoscopic Instrument Control Systems (da Vinci X and da Vinci Xi Surgical Systems) are intended to assist in the accurate control of Intuitive Surgical Endoscopic Instruments during urologic surgical procedures, general laparoscopic surgical procedures, gynecologic laparoscopic surgical procedures, general thoracoscopic surgical procedures, and trans-oral otolaryngology surgical procedures restricted to benign tumors and malignant tumors classified as T1 and T2, and for benign base of tongue resection procedures. The systems are indicated for adult and pediatric use (except for trans-oral otolaryngology surgical procedures). They are intended to be used by trained physicians in an operating room environment.

The da Vinci X and da Vinci Xi Surgical Systems are class IIb medical devices CE marked (CE 2460) under the European Medical Devices Directive (93/42/EEC), manufactured by Intuitive Surgical, Inc. Refer to Instructions For Use before use.

SureForm 45 & 60 Staplers

The SureForm 45 Stapler & SureForm 60 Stapler and SureForm 45 Stapler Reloads & SureForm 60 Stapler Reloads and other Stapler Accessories (including the bladeless obturators) are intended to be used with a compatible da Vinci Surgical System for resection, transection, and, or creation of anastomoses in General, Thoracic, Gynecologic, Urologic, and Pediatric surgery. The devices can be used with staple line or tissue buttressing material (natural or synthetic).

The SureForm 60 Stapler and SureForm 45 Stapler, SureForm 60 and SureForm 45 Stapler Reloads are class IIa and IIb medical devices CE marked (CE 2460) under the European Medical Devices Directive (93/42/EEC), manufactured by Intuitive Surgical, Inc. Refer to Instructions For Use before use.

8 mm SureForm 30 / SureForm 30 Curved-Tip

The Intuitive Surgical 8 mm SureForm 30 Staplers, 8 mm SureForm 30 Reloads, and accessories are intended to be used with a compatible da Vinci Surgical System for resection, transection of vasculature and tissue, and/or creation of anastomoses in General, Thoracic, Gynecologic, Urologic, and Pediatric Surgery.

The Intuitive Surgical 8 mm SureForm 30 Staplers, 8 mm SureForm 30 Reloads, and accessories are class IIa and IIb medical devices CE marked (CE 2460) under the European Medical Devices Regulation (EU 2017/745), manufactured by Intuitive Surgical, Inc. Refer to Instructions For Use before use.

Intuitive's Privacy Notice is available at www.intuitive.com/privacy.

© 2024 Intuitive Surgical Operations, Inc. All rights reserved. Product and brand names/logos are trademarks or registered trademarks of Intuitive Surgical or their respective owner.

