

# Da Vinci Handheld Camera

Begin building an efficient operating room by simplifying setup. Perform da Vinci cases from start to finish without relying on a separate laparoscopy tower. Accomplish optical first entry through an 8 mm port.



**Compatible with all third-party 2D laparoscopes with a standard eyepiece (ISO 18339:2015)\***



\*Any third-party 5–10 mm rigid laparoscope with an eyepiece cap compliant with ISO/TS 18339 can be connected to the da Vinci Handheld Camera.

## Ordering Information

| Part number | Product Description                 | Compatibility                          |
|-------------|-------------------------------------|--|
| 470035      | Da Vinci Handheld Camera            | da Vinci X, Xi and SP Surgical Systems |
| 470042      | Handheld Camera Light Guide         | da Vinci X, Xi and SP Surgical Systems |
| 470051      | Handheld Camera Light Guide Adapter | da Vinci X, Xi and SP Surgical Systems |
| 400499      | Handheld Camera Sterilization Tray  | da Vinci X, Xi and SP Surgical Systems |

**To learn more about da Vinci Handheld Camera, visit:**  
[www.intuitive.com](http://www.intuitive.com)

**To contact a representative or receive additional information:**

Visit [www.intuitive.com](http://www.intuitive.com) or call Intuitive Customer Service: +41.21.821.2020 or +800.0821.2020 (toll free)

Medical device, CE, refer to Instructions For Use for further information.

Individuals' outcomes may depend on a number of factors, including but not limited to patient characteristics, disease characteristics and/or physician/surgeon experience.

Some products, features or technologies may not be available in all countries. Please contact your local Intuitive representative for product availability in your region. Refer to

the product specific User Manual for indications, contraindications, warnings and other product information.

Intuitive's Privacy Notice is available at [www.intuitive.com/privacy](http://www.intuitive.com/privacy).

© 2020-2024 Intuitive Surgical Operations, Inc. All rights reserved. Product and brand names/logos are trademarks or registered trademarks of Intuitive Surgical or their respective owner.