

Da Vinci X

Expand your
potential.



INTUITIVE



Expand your potential

As a cost-effective gateway to the fourth generation of da Vinci surgical systems, the da Vinci X can grow with your robotics program as you perform innovative minimally invasive surgery for your patients.

Designed so you can comfortably perform a range of procedures, this system provides visual guidance during setup, a linear port placement, and a smaller, thinner arm design. This enables a broader reach than previous generations of da Vinci systems.¹

Plus, you can use much of the same advanced technology available for our most capable platform, the da Vinci Xi surgical system, making upgrades straightforward if your needs evolve.

Build on a solid foundation

Intuitive created the da Vinci X surgical system by using experience from three previous generations and thousands of hours working with surgeons and staff. Now the benefits of robotic-assisted surgery can be brought to any operating room at a lower price.²

Set up quickly

Follow a clear interface as you set up, and utilize a laser guidance system that helps take the guesswork out of positioning the patient cart.

Simplify port placement


The architecture of the fourth-generation instrument arms allows ports to be placed in a straight line for multi-port procedures.

Switch your view

A universal arm design allows you to insert the endoscope on any arm at any time during your procedure.


Get a greater range of motion¹

By combining longer instruments with thinner arms¹, the da Vinci X surgical system can help limit the potential for collision while focusing on precision.



Access the latest da Vinci innovations

Choosing the da Vinci X surgical system means you can use the same instruments available for the da Vinci Xi, enabling you to standardize instruments for your multi-port patient carts.



See tissue in a different light

Built-in Firefly fluorescence imaging gives you real-time visualization and assessment of vessels, bile ducts, and tissue perfusion.



Put the power in your hands

The integrated ERBE VIO dV generator helps you fine-tune energy output and switch presets when you need to directly from the surgeon console.



Staple with control

Da Vinci Stapling offers a full range of lengths and reloads—30, 45, and 60 mm—as well as fully wristed articulation, complete surgeon control from the console, and software-driven feedback to help reduce guesswork.



Get two grip strengths in one grasper

Force Bipolar with DualGrip technology can help streamline setup and minimize workflow disruptions by reducing instrument exchanges.



MEDLINE

Upgrade when you are ready

Powered by da Vinci OS4, the fourth generation of da Vinci surgical systems allows flexible configurations and an upgrade path as your program needs evolve, whether you are switching to the da Vinci Xi or adding future technologies.



System comparison

Da Vinci X



Da Vinci Xi



Core technology

- Da Vinci OS4 compatible
- Intuitive motion
- EndoWrist technology
- Immersive 3DHD stereo viewer
- Upgradeable architecture



Technology enhancements¹

- Dual console compatible
- Single-Site technology
- Vessel Sealer Extend
- Firefly Fluorescence Imaging
- Multi-Port instruments and endoscope
- Integrated energy
- Da Vinci stapling portfolio
- Platform modularity
- Handheld camera
- Endoscope Plus



Capability enhancements¹

- Streamlined port placement
- Extended anatomical access
- Integrated Table Motion
- Boom rotation



Ease of use enhancements¹

- Automated setup tasks
- Advanced intraoperative adjustments
- Universal cart positioning
- Optimized patient side access



Purchasing information

Da Vinci X system

To contact a representative or receive additional information visit:
www.intuitive.com

or

call Intuitive Customer Service in India at:

+180 0103 6952 (toll-free) / +91 80 455 05 100




What you get³

The da Vinci X system has three components:

1. A unique **patient cart**

2. A **surgeon console** providing an immersive, ergonomic experience

3. A **vision cart** supporting the latest 3DHD vision system



Important safety information

Serious complications may occur in any surgery, using da Vinci surgical system, up to and including death. Examples of serious or life-threatening complications, which may require prolonged and/or unexpected hospitalization and/or reoperation, include but are not limited to, one or more of the following: injury to tissues/organs, bleeding, infection and internal scarring that can cause long-lasting dysfunction/pain.

Risks specific to minimally invasive surgery, using da Vinci surgical system, include but are not limited to, one or more of the following: temporary pain/nerve injury associated with positioning; a longer operative time, the need to convert to an open approach, or the need for additional or larger incision sites. Converting the procedure could result in a longer operative time, a longer time under anesthesia, and could lead to increased complications. Contraindications applicable to the use of conventional endoscopic instruments also apply to the use of all da Vinci instruments.

For Important Safety Information, indications for use, risks, full cautions and warnings, please also refer to www.intuitive.com/safety.

Individuals' outcomes may depend on a number of factors, including but not limited to patient characteristics, disease characteristics and/or surgeon experience.

Da Vinci Surgical Systems

The Intuitive Surgical Endoscopic Instrument Control Systems (da Vinci X, da Vinci Xi and da Vinci Si surgical systems) are intended to assist in the accurate control of Intuitive Surgical Endoscopic Instruments during urologic surgical procedures, general laparoscopic surgical procedures, gynecologic laparoscopic surgical procedures, general thoracoscopic surgical procedures, thoracoscopically-assisted cardiomy procedures, and trans-oral otolaryngology surgical procedures restricted to benign tumors and malignant tumors classified as T1 and T2, and for benign base of tongue resection procedures. The systems can also be employed with adjunctive mediastinotomy to perform coronary anastomosis during cardiac revascularization. The systems are indicated for adult and pediatric use (except for trans-oral otolaryngology surgical procedures). They are intended to be used by trained physicians in an operating room environment.

The da Vinci X, da Vinci Xi and da Vinci Si surgical systems are class IIb medical devices CE marked

(CE0543) under the European Medical Devices Directive (93/42/EEC), manufactured by Intuitive Surgical, Inc. Refer to Instructions For Use before use.

Firefly Fluorescence Imaging

The da Vinci Firefly imaging system is intended to provide real-time endoscopic visible and near-infrared fluorescence imaging. The da Vinci Firefly imaging system enables surgeons to perform minimally invasive surgery using standard endoscopic visible light as well as visual assessment of vessels, blood flow and related tissue perfusion, and at least one of the major extra-hepatic bile ducts (cystic duct, common bile duct or common hepatic duct), using near infrared imaging.

Table motion

The da Vinci table motion is intended to allow the surgical staff to reposition the patient by adjusting the table without undocking the da Vinci Xi surgical system during urologic surgical procedures, general laparoscopic surgical procedures, and gynecologic surgical procedures. It is designed to be used with a compatible operating room (OR) table.

SureForm staplers

The SureForm 60 and SureForm 45, SureForm 60 and SureForm 45 reloads are intended to be used with a compatible da Vinci surgical system for resection, transection, and, or creation of anastomoses in general, thoracic, gynecologic, urologic, and pediatric surgery. The devices can be used with staple line or tissue buttressing material (natural or synthetic).

The SureForm 60 and SureForm 45, SureForm 60 and SureForm 45 reloads are class IIa and IIb medical devices CE marked (CE0543) under the European Medical Devices Directive (93/42/EEC), manufactured by Intuitive Surgical, Inc. Refer to Instructions For Use before use.

Force Bipolar

The EndoWrist Force Bipolar is intended to be used with compatible systems for endoscopic manipulation of tissue, including dissection, grasping, retraction, and bipolar coagulation of tissue.

The EndoWrist Force Bipolar is a class IIb medical device CE marked (CE0543) under the European Medical Devices Directive (93/42/EEC) manufactured by Intuitive Surgical, Inc. Refer to Instructions For Use before use.

Vessel Sealer Extend

The Vessel Sealer Extend is a bipolar electro-surgical sealing and cutting instruments for use with the da Vinci Xi and da Vinci X surgical systems and the ERBE VIO dV electro-surgical generator. It is intended for grasping and blunt dissection of tissue and for bipolar coagulation and mechanical transection of vessels up to 7 mm in diameter and tissue bundles that fit in the jaws of the instrument. The Vessel Sealer Extend has not been shown to be effective for tubal sterilization or tubal coagulation for sterilization procedures, and should not be used for these procedures.

The Vessel Sealer Extend is class IIb medical devices CE marked (CE0543) under the European Medical Devices Directive (93/42/EEC), manufactured by Intuitive Surgical, Inc. Refer to Instructions For Use before use.

ERBE VIO dV

The ERBE VIO dV with accessories electro-surgical unit is intended to deliver High Frequency (HF) current for the cutting and/or coagulation of tissue.

The ERBE VIO dV is a class IIb medical device CE marked (CE 0124) under the European Medical Devices Directive (93/42/EEC) manufactured by ERBE Elektromedizin GmbH. Refer to Instructions For Use before use.

Unless otherwise noted, all people depicted are models.

References

1. Compared to the da Vinci Si surgical system. Results of testing on file.
2. Lower price refers to initial expenditure on surgical system.
3. Parts of the da Vinci X surgical system may be refurbished or reconditioned.

PN 1114574_RevA_IN_02/2023

*For registered medical practitioners use only

India Headquarters & Registered Office

Intuitive Surgical India Private Limited

The Millennia, Tower A, Level 5,
No. 1 & 2 Murphy Road, Ulsoor
Bengaluru – 560 008

CIN: U74120KA2014FTC146332

Tel: + 91 80 455 05 100
Fax: +91 80 455 05 126

Global headquarters

USA

Intuitive Surgical, Inc.
1020 Kifer Road, Building 108
Sunnyvale, CA 94086-5304
Tel: +1.408.523.2100

For more information, call Intuitive

+180 0103 6952 (toll-free) / +9180 455 05 100
or visit www.intuitive.com

INTUITIVE

